

SOP 19 *In Vitro* Cell Culture Expansion Procedure

Preparation

- 1) Wash hands, glove up; aseptically clean the bio-safety cabinet (hood) with 70% ethanol and sterile wipes
- 2) Warm medium to 37^o C in the water bath
- 3) Clean 1 bottle of PUCKS EDTA with 70% ethanol and place into hood
- 4) Gather flask(s) to be expanded and place into hood
- 5) Label new centrifuge tube(s) with cell line name
- 6) Using a pipette, remove all media+cells from the flask(s) and pipette into corresponding labeled centrifuge tubes
- 7) Grab a new pipette and add PUCKS into the flask(s) with a "sweeping motion" to evenly distribute the solution into the flask.
 - a. Place flask(s) back into the correct incubator for 10-15minutes to allow time for the cells to detach
- 8) While waiting:
 - a. Get the media from the hot bath; clean with 70% ethanol and place into hood
 - b. Clean area
- 9) After 15 minutes, get the flask(s) from the incubator, tap the side of the flask with your hand to loosen cells and place into the hood
- 10) Pipette the PUCKS+cells into the same labeled centrifuge tubes used in step 6.
- 11) Centrifuge all but one tube (keep one set aside for combining and counting later in step 13) at 500 x g (1600rpm) for 6 minutes
- 12) After 6 minutes, pour off supernatant as waste without disrupting the new cell pellet
- 13) Using the set aside (unspun) centrifuge tube with media+cells, combine all cell pellets. Mix well
- 14) Proceed to count the cells

Counting Cells using Vi-Cell XR

- 1) Get a Vi-Cell™XR counting vial and add 500 µL of the mixed sample
- 2) Record the total volume of the centrifuge tube
- 3) Count the cells using the Vi-Cell™XR
 - a. Click: Log in sample
 - b. Enter: Sample ID, Cell Type, Dilution Factor (1)
 - c. Click: Ok and start queue
- 4) Multiply total viable cells by the total volume and record the cell viability
- 5) Calculate how many cells are needed for each experiment moving forward (ie. those cells going back to culture, being pelleted, cryopreserved and/or for other purposes – each will need a separate calculation to determine cells/volume needed)
 - a. *Calculation: (Total cells needed ÷ Total viable cells) x total volume = volume needed to aliquot cell count needed*
- 6) Once all calculations are complete, label new centrifuge tubes with cell line name and purpose of the tube.
- 7) Aliquot appropriate media+cell volumes to distribute the cells accordingly
- 8) Centrifuge all tubes at 500 x g (1600rpm) for 6 minutes, pour off supernatant as waste
- 9) Proceed to each appropriate next step

Cryopreservation

- 1) Cryopreservation materials:
 - a. Outer-threaded vials
 - b. 1 tube of pure FBS
 - c. 1 tube of 15% DMSO in RPMI-1640 medium
 - d. Isopropanol tub

- e. Centrifuge tube containing cells for freezing
- 2) For every vial cryopreserved, add 0.8mL (800 μ L) of FBS and 0.8 mL (800 μ L) of 15% DMSO in RPMI-1640 medium directly to the centrifuge tube containing the cells to be frozen.
 - a. Final concentrations of cryopreservation reagents are: 50% FBS and 7.5% DMSO
- 3) Mix with the serological pipette
- 4) Add 1.6 mL of the viable cells suspended in the cryopreservation medium + FBS to each 2 mL outer-threaded vial being frozen
 - a. Label vials with: Cell line name, passage, cell count, viability, "Geno/MycO OK" today's date and your initials.
- 5) Place vials into the isopropanol control rate freezing tub and into a -80°C freezer
 - b. Samples are to be inventoried and transferred 12-24 hours later into a liquid nitrogen vapor freezer
- 6) Clean area with 70% ethanol

Pelleting

- 1) Re-suspend the cell pellet in PBS to wash the cells
- 2) Centrifuge at 500 x g (1600rpm) for 6 minutes, pour off PBS supernatant as waste
- 3) Re-suspend the cell pellet in 0.5mL PBS per vial
 - a. If you are pelleting 5 vials at 5e⁶ cells each, re-suspend the 25e⁶ cell pellet in 2.5mL PBS total
- 4) Transfer 0.5mL of PBS+cells into each labeled 1.5mL microcentrifuge tube.
 - a. Label vials with: Cell line name, passage, cell count, viability, "Geno/MycO OK" today's date and your initials.
- 5) Using a microcentrifuge, centrifuge the vials at roughly 3000-4000rpm for 1-2 minutes
- 6) Pipette off the PBS supernatant as waste
- 7) "Snap freeze" the dry pellets in liquid nitrogen for 10 seconds each and immediately store at -80°C
- 8) Clean area with 70% ethanol

Re-Culturing Cells

- 1) Obtain desired number of flask(s)
- 2) Label flask(s): Cell line name, passage (increase by 1), today's date, and your initials
 - a. Also add any additional notes necessary to your experiment to your flask(s)
- 3) Re-suspend (mix well) the cell pellet in the appropriate amount of media (see below) and place into the base of each newly labeled flask
 - a. Amount of media will vary according to the number and size of the flasks you're expanding into
 - b. Flask sizes and media amounts: T12.5 = 5mL, T25 = 7mL, T75 = 13mL, T150 = 25mL, 3-layer = 96mL, 5-layer = 160mL
- 4) Place the flask(s) back into their corresponding incubators
- 5) Clean area with 70% ethanol

Maintaining Cultures:

- 1) At the start and end of each experiment, submit a 2e6 dry cell pellet for STR analysis and 1mL of supernatant (no cells) for mycoplasma assay analysis.
- 2) Flasks need to be routinely checked via a microscope every 7-10 days.
 - a. Incubators need to be monitored daily to verify that no flask is visibly contaminated or if the media needs replenishing earlier than 7-10 days.
 - b. Watch for cell proliferation, contamination, and growth of non-cancerous cells via the microscope.
- 3) Medium needs to be changed every 7-10 days to maintain fresh nutrient concentrations
- 4) Flask(s) should be expanded when 80-90% confluent or when optimal number of cells needed for a specific experiment are apparent
- 5) Discard the cells immediately when microbial contamination is present.